

Ermes-Tipco Order form



Order Quotation

Invoicing Details	Shipping Details				
Company name:	Ship to:				
Delivery date required:/					
Delivery address:	CUSTOMER STAMP AND SIGNATURE				
P. Code.: Town:					
Tel.: Fax:					
# Order: Date of the Order:/					

Pro Medicare srl grants **24 months** warranty on production defects from the date of the first put into service and **12 months** on covers and components subjected to wearing. Warranty is void if the product is used in any other way than stated in the Instruction Manual.


Note: Any request for variation in dimensions, materials and type compared to the standard does not have **CE** marking, therefore the professional user becomes the manufacturer and has the obligation to draw up the technical documentation of the device. *Images are purely illustrative and may not fully reflect reality.*

Ermes-Tipco Standard Configuration

- Seat adjustable in width, depth also asymmetrically in adduction and abduction
- Kit of inserts for a customised positioning of the pelvis
- Reclining and height adjustable backrest
- Kit of inserts for a customised positioning of the trunk
- Multi-adjustable and removable containment and positioning kit:
 - pelvis (adjustable in width, depth, height, adduction/abduction)
 - arm support (adjustable in width, depth, height, angle, intra/extra rotation)
- Male interface plate
- Tipco indoor frame adapted with female interface plate:
 - Height adjustable seat by gas spring
 - Tilt in space by gas spring from -5° up to ° +30°
 - 4 swivel castors (2 with locks)
 - Multi-adjustable push handle



Sizes of the user



1 RH	depth of the seat RH	
1 LH	depth of the seat LH	
2 RH	length of the leg RH	
2 LH	length of the leg LH	
3 RH	depth of the thorax RH	
3 LH	depth of the thorax LH	
4	height of the shoulders from the seat plane	
5	width of the thorax	
6	width of the pelvis	

Ermes

	•					
0	size 2	color	0	blue/grey	0	fuchsia/grey
0	size 3	color	0	blue/grey	0	fuchsia/grey
0	size 4	color	0	wine/grey		

Note: If Ermes' sizes 2/3 are choosen, Tipco indoor base will be in white color; if Ermes' size 4 is choosen, Tipco indoor base will be in grey color

Note: $O = optional without up-charge \square = optional with up-charge$

Mod.Ord.Erm.Tip. S.P. EN Rev.7 10/2023 1/4

Type of the seat

MIDI Seat

lateral/anterior control (standard configuration)



DEEP Seat total control (PTS) <u>optional without</u> <u>up-charge</u>



				100	District Control of the Control of t							and the second				
□ Grea	ter closure	of the se	at-ba	ckrest a	angle	(±5°)										
Thorac	ic Suppor	<u>ts</u>														
□ Flat	pads								□ Con	toured pac	 Is					
7		Size	(cm)	1	2	3	4			A	Size (cm)	1	2	3	4	
H		Н		7	7	8	10		H		Н	6	7	8	10	
~	s= 2,5 cm	L		7.5	8.5	10	12		~	s= 2,5 cm	L	9.5	11	13	16	
Side	Size	Bra	cket	type*		Bracket shape			Side	Size	Bracket	Bracket type*			Bracket shape	
□ RH		□ Fixed					Α		□ RH		□ Fixed				A	
		☐ Swing-A	way				В				☐ Swing-Away				В	
□ LH		□ Fixed					Α		□ LH		□ Fixed				A	
		□ Swing-A									☐ Swing-Away					
	echnical infori ffset Bracket)		e end o	of this for	m you	can find th	e dimensi	ons	reachable	with the two	types of brackets (A	4 type = l	_/90° Bra	cket; B ty	/pe =	
Belts a	nd Harne	sses (value	s expre	essed in c	m, see	details in t	he techni	cal	informatio	n)						
Ches	t pad					Dynamic	Adbuc	tio	n Harnes	ss	□ Pelvic belt					
	00	Size		HxW			-			XS (2)*			Size	н	xW	
8.	9	□ 1	(8x10)			2			S (3)*	PERM	SI	□ 1	(6:	x11)	
		□ 2	(9x13)		138		2		M (4)*			□ 2	(7:	x13)	
	racic supports only)	□ 3	(1	L0x13)		V.							□ 3	(7	x15)	
Dyna	mic Butter	fly Harnes	ss		_ J	Jacket Ha	arness				□ Contoured	Pelvic	Belt			
-	,		х	S (2)*		1	4			XS (2)*			Size	Н	xW	
			9	S (3)*			ì			S (3)*	100	1	□ 1	(6:	x11)	
			N	/ 1 (4)*		H				M (4)*	00		□ 2	(8:	x15)	
														I is referre of the cent of the pad	tral part	
*The nui	mber in the b	rackets is ref	erred	to the rele	evant E	rmes size.	Therefore	e it i	is advisabl	e to combine	the default size.	l				
Hardre	est with m	ulti-axial j	oint a	nd muli	ti-adju	ustable a	rm									
Confort						atomic v		1			3 piece pads			eadrest	with	
	chable hard hardware	aware				detachal fixed har		lWa	are		☐ detachable		vare			
				S	1					XS			□ 1	(7:	x13)	
				M						S			□ 2	(9:	x20)	
				G	•					M				expressed i		
Other	accessoi	ries													•	
			(n	nee elev	for an <u>c</u>	se	parate fo	oot	□ LH plates:		H					

Mod.Ord.Erm.Tip. S.P. EN Rev.7 10/2023 2/4

☐ footboard (it does not allow asymmetry, adduction and abduction)



□ Shoes Holders



Calf support □ RH □ LH

□ **Calf panel** (it does not allow asymmetry, adduction and abduction)


Lateral supporting pad for legs □ RH □ LH


- □ Padded anterior/posterior Calf strap
- □ Padding for separate footplates
- □ Padding for footboard


Elbow lateral-posterior support □ RH □ LH

- ☐ Grey detachable Tray with body contouring recess and edges
 - O standard version
 - O larger version for daily activities
- ☐ Transparent detachable Tray with body contouring recess and edges (only standard version)

Note: the dimensions of the tray in the <u>standard</u> grey color version and those of the transparent tray are the same, the naming changes (the size 2 of the grey tray will be the size 1 of the transparent tray, 3->2, 4->3)







Size (cm)	L	P	P1
1 (for Ermes size 2)	50	45	13,7
2 (for Ermes size 3)	60	50	15.4
3 (for Ermes size 4)	68	55.6	19.3

Padding of the tray □ upper

□ lower

□ Addictional Outdoor Frame Adacta Giò for Ermes

Standard Configuration

- Adacta titanium satin anodised frame which plates have the same color of the cover
- Adjustable single handlebar with folding backrest tubes
- Reclining mechanism by predetermined plate w/pin 3°, 10°, 18°
- Continuous Tilt in space by drive lever: from 2° up to 35° (from 6° up to 39°) tilt in space can be up-graded
 of 4° with an increment of the starting inclination of the frame by operating on the positioning holes of the casters forks
- Fixed anti-tip system
- Wheels Kit with brake
- Height adjustment: high (49 cm from the floor)/low (43 cm from the floor) with 300mm and 500mm wheels only
- Solid grey base assembled with female interface plate



◊ = Standard	= Optional without	t up-charge □ = 0	Optional with up-charge	2			
Frame size	2 (combined w	vith Ermes size 2)	3 (combined	with Ermes size 3)	6 (combined with Ermes size 4)		
<u>Castors</u>	Castors ♦ 150mm (6") poly grey ○ black ○		O 175mm (7") poly	O grey	O black		
Ruote posteriori	O 300mm (12") p O 300mm (12") p O 300mm (12") p	pooly grey with park poly black with park poly black with drund poly grey with drund poly black with drund poly grey with drund poly grey with drund poly black with drund poly black with park poly black with drund poly black with drund poly grey with drund poly black with drund poly black with drund poly grey without he po	king brake lever king brake lever h parking brake lever m brake h drum brake m brake and parking br h drum brake and parking br h drum brake and park king brake lever (only m brake (only 45 cm) m brake and parking br andrims with parking b	ake lever ake lever ing brake lever 45 cm height) height) ake lever (only 45 cm in the second control of the s		O BIUCK	
	□ 500mm (20")	pneumatic <i>grey</i> wi	th handrims with drum	brake and parking bral	ke lever		

Mod.Ord.Erm.Tip. S.P. EN Rev.7 10/2023 3/4

	O Kit for fixed push handles								
	☐ Kit for telescopic and angle adjustable single handlebar								
Back Posts	☐ Kit for telescopic push handles with lever block (max adjustment 10 cm)								
	☐ Kit for angle adjustable separated handlebar								
	☐ Quick-release pins for easier back reclining								
Tilt in space	vith pedal O with lever								
Canopy adaptor in	cluded Blue Pink Bordeaux								
□ <u>Basket</u>									
□ Rain cover									
Foot muff	□ Blue □ Pink □ Bordeaux								

TECHNICAL INFORMATION

Ermes		Size 2	Size 3	Size 4	
Width of the Seat (cm)	Width of the Seat (cm)				
Depth of the Seat (cm)	Depth of the Seat (cm)				
Height of the trunk (cm)	26-38	33-45	41-53		
\Midth of the therew (am)	bracket A 19-25	bracket A 22-28	bracket A 27-33		
Width of the thorax (cm)		bracket B 14-20	bracket B 17-23	bracket B 22-28	
Max abduction of the leg (°)		16	16	16	
Max adduction of the leg (°)		16	16	16	
Popliteal Space to footplate Length (cm)* *if the elevating knee angle is provided		8-19	8-29	8-39	
Height of the arm support (cm)		14-19	20-26	20-26	
Backrest Reclining (°)		-3+20	-3+20	-3+20	
Ermes load (kg)		25	35	50	
Ermes weight (kg)		5.2	8.1	12.1	
Contoured hip guide		Size 2	Size 3	Size 4	
н Д	Ι				
	Н	6.5	7	8	
p	P	17	20	23	
Flat arm support		Size 2	Size 3	Size 4	
	L	7.5	7.5	8.5	
P	Р	17	20	23	
Note: the postural unit for armrest and the lateral hip pad are	comb	ined to the relative size o	of Ermes and are not inte	rchangeable	
Tipco		Size 2	Size 2	Size 3	
Seat to floor height (cm)		63	63	63	
Width overall dimension (cm)		67	67	67	
Depth overall dimensions (cm)		80	80	88	
Tilt in space (°)		-5+30	-5+30	-5+30	
Tipco Load (kg)		50	50	50	
Tipco Weight (kg)		13.6	13.6	14	
Total Weight (<i>Tipco+Ermes</i>) (kg)		18.8	21.7	26.1	
Adacta Giò		Size 2	Size 3	Size 6	
Seat to floor height (cm)		43	43	43	
Width overall dimension (400/500 wheels) (cm)		48	53	58	
Min/max depth overall dimensions (cm)	73-100	75-102	81-109		
Tilt in space (°)	2-35 (from 6 up to 39) tilt in space can be up-graded of 4° with an increment of the starting inclination of the frame by operating on the positioning holes of the casters forks				
Predetermined backposts reclining (°)		3-10-18			
Adacta Giò load (standard configuration) (kg)		30	50	50	
Adacta Giò weight (kg)		12.8	13.4	14.5	
Total weight (Adacta Giò+Ermes) (kg)		18	21.5	26.6	



PRO MEDICARE S.r.l.

Via Montagna Z.I. Lotto 41-72023 Mesagne (Br)-ITALY Tel: +39-0831-777840 Fax:+39-0831-730739 Website: www.promedicare.eu e-mail: sales@promedicare.it

Mod.Ord.Erm.Tip. S.P. EN Rev.7 10/2023 4/4